

LA049179UFEUB7

Antibac Telescopic Mop And Refill

Please retain instructions for future reference.

CARE AND USE

To be eligible for the guarantee, please retain your proof of purchase.

Please retain instructions for future reference.

Directions for use Before first use, remove all of the packaging.

Extend the mop handle to a comfortable height and twist it in a clockwise direction to lock it into position.

Replacing the mop head Twist the mop head in an anticlockwise direction to unlock it.

Remove the used mop head and then replace it by pushing the refill mop head onto the end of the mop and twisting in a clockwise direction until it locks into place.

Caution Use only the recommended amount of detergent or disinfectant, according to the manufacturer's instructions.

Dimensions stated are approximate.

Made in China

Manufactured by: UP Global Sourcing UK Ltd, Victoria Street, Oldham, Manchester OL9 0DD. UK
Ultimate Products Europe Ltd, 19 Baggot Street Lower, Dublin, D02 X658. ROI